


SITING REQUIREMENTS



- ▶ Sites within a minimum of 550 sq. ft.
- ▶ Weighs approx. 8,500 lbs.
- ▶ Short bore active shielded super-conducting magnet
- ▶ R/F Shielding: > 100 db attenuation
- ▶ Power: Main service 120/208/150A/3 Phase/60 Hz
- ▶ Air Conditioning (Magnet and Equipment Rooms): 9,200 BTU/hr. min.
- ▶ Temperature range: 68-74°F in scan room
- ▶ Humidity range: 30-60% relative